



Recombinant Mouse PWWP domain-containing protein 2A (Pwwp2a)

Product Code	CSB-BP731676MO
Abbreviation	Pwwp2a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q69Z61
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MAAVAAEAAA TAASPGEGGA GEAPELEPI PGSEAGTPLP VTATEAAVPD GEADGRQSAP QADEQPLPPP PPPPPPGELA DSSEAEAAKP PEPAAVPVSP PEQPPAAPEQ PEDAPRPPPA PALVPPAGGD SAVSHLIPGS EVRVTLDHII EDALVVSFRL GEKLFSGVLM DLSKRFGPHG IPVTVFPKRE YKDKPDAMQL QSTTFQEGIE VKQEVNGAVP DDLSPVPPPE RLWASKPPPL FHEGAPYPPP LFIRDTYNQS IPQPPRRIK RPKRKMREE PTSIMNAIKL RPRQVLCDC KNSVVAEKKE IRKGSSDSSR YEDKKRRNDS VATVNKKLKT DHKVDGKNQN ESQRRNTVVR VSSIAHSRGR VVKVSAQANT SKAQLNKKV LQSKNMDHAK AREVLKIAKE KAQKKQSETS TSKTAHAKVH FTRRYQSPSS GSLPPRVRLK PQRYRNEEND SSLKTGLEKI RSGKLAPKPQ SRCTSTRSAG EAPSEKPSPS EGPEESAGEV QDTSRVRVPG EQEELRMLGK KGSKSSISVY LTLNQETSDS SSASVCSIDS MDDLKSSNSE CSSSEFVFP PGCMHAPSAS STSTSFSSKE ENSLRNSLKM KIFSKNVSKC ITPDGRTICV GDIVWAKIYG FPWWPARILT ITVSRKDNGL LARQEARISW FASPTTSSLA LSQSPFLEN FQLRFNKKRK GLYRRAITEA AKAQKLTPE VRALLTQFET</p>
Source	Baculovirus
Target Names	Pwwp2a
Protein Names	Recommended name: PWWP domain-containing protein 2A
Expression Region	1-730
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.