



Recombinant Mouse Raftlin (Rftn1)

Product Code	CSB-EP731683MO
Abbreviation	Rftn1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6A0D4
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	<p>GCSLNKLEK REEKRPGNIY STLKRPQVET KVDVTYEYCF LEFTTLTAAE LPRSSATRLA SLRDLPDQLL ELYQQGFSLA ALHPFVQPTR RQEKILLEHI FRAILVKKTN RSQKAELHDE GYTLLELDYCS SLEHLADQKL IPEFIKKVQE AASQGLKFVS VVPQYQPSVS SAGSRRLKPV ANSVEDARDV KCTLGDRSSL ENDTPKAAET DAATAGVNRN PETKPGSTGD VPSAQQPGIP SPSAENEAGE FPLRGLQPAL DRSEGDPNSG PEELPSRKME IFAFFNRPKS QQKCRQYYPV TIPLQVSKNG QTVSSLDASW LEHMSDHFRK GGVLVNAV FQ LGMANDSFYQ LTDGVFIFEAVSTEDNRTTQ GYDAIVVEQW TVLEGTEVQT DYMPLLNSLA AYGWQLTCVL PTPILKTTRE GNVSTKQIVF LQRPCLPQKT KKRESKFQWR FSRNEIHGRQ TRKSKGKLSA SNKQQAEEENE KNLEDQFSKA GDVGNLVLGA PQWGRASEVR EQRQGSAAVQ NGPAGHNRDS VALRHSNPRA EAELAAGPTP TEAN</p>
Source	E.coli
Target Names	Rftn1
Protein Names	Recommended name: Raftlin Alternative name(s): Raft-linking protein
Expression Region	2-554
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.