



Recombinant Propionibacterium acnes Peptide chain release factor 2 (prfB)

Product Code	CSB-MP731704PSG
Abbreviation	prfB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6A808
Product Type	Recombinant Protein
Immunogen Species	Propionibacterium acnes (strain KPA171202 / DSM 16379)
Purity	≥85% (SDS-PAGE)
Sequence	MALVDLQQTV EDLDASLSSI ENVLDPEAKK AEIADLERQV AVPNLWDDQE NAQRVTSRLS ALQAEIERLE RLRARIEDVH VLLEFAAADD DKESLAEATT ETVKLREEID ALEVRTLSSG QYDEREAVVT IRSEAGGVDA ADFAEMLLRM YLRWSERHGY KVEILDTSYA EEAGIKSATF VVHAPFAYGT LSVEQGTHRL VRISPFDNQG RRQTSFAGVE VLPVVEEDH VDIPESDLRV DVFHASGPGG QGVNTTDSAV RLTHLPTGIV VSCQNERSQI QNKAAALRVL QAKLLEKARE EQEAEMNSLK SEGNSWGAQM RSYVLHPYQM VKDLRTNHEV GNTDAVFDGD IDDFIDAGIR WRRKNETD
Source	Mammalian cell
Target Names	prfB
Protein Names	Recommended name: Peptide chain release factor 2 Short name= RF-2
Expression Region	1-368
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.