



Recombinant Leifsonia xyli subsp. xyli Phosphoribosylaminoimidazole- succinocarboxamide synthase (purC)

Product Code	CSB-YP731725LAAC
Abbreviation	purC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ACE5
Product Type	Recombinant Protein
Immunogen Species	Leifsonia xyli subsp. xyli (strain CTCB07)
Purity	≥85% (SDS-PAGE)
Sequence	MTEMPGWKHV YSGKVRDLFV PETAEGLTET DTVLVVASDR VSAFDHVLEP GIPGKGELLT ALSLWWFDRL DAPNHLIPDH TLEGERARER IPAAVSGRSM LVKPLDMFPI ECVVRGYLTG SGWQEYQESQ SVCGLPLPAG LSDGDRLPEP IYTPAWKAPL GEHDENISYE RTIELVGAET AGTLRELSLA VYGRGAAIAE KRGVIIADTK FEFGAERETG EITLADEVLT SDSSRYWDAA ATSDVYTPKE RMASFDKQIV RDWLAANWDR TGTPPALPAE IVERTAGRYR ALVERLTGA
Source	Yeast
Target Names	purC
Protein Names	Recommended name: Phosphoribosylaminoimidazole-succinocarboxamide synthase EC= 6.3.2.6 Alternative name(s): SAICAR synthetase
Expression Region	1-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.