



Recombinant *Desulfotalea psychrophila* NADH-quinone oxidoreductase subunit D (nuoD)

Product Code	CSB-BP731779DAAA
Abbreviation	nuoD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ANM7
Product Type	Recombinant Protein
Immunogen Species	<i>Desulfotalea psychrophila</i> (strain LSv54 / DSM 12343)
Purity	≥85% (SDS-PAGE)
Sequence	MSEHIILNPD ETFTLNLGPQ HPATHGVLRV KLTMDGEYIY SAETVIGYIH RMHEKMGENR TYNQYLPNLS RLDYLSAMAY CHGYVLAVEK AGSIEVPERA EYIRAITVEL NRLSSHLVWF GAFIMDLGGF SPLMYAFDDR EQILDLLESI TGSRLTYCYY RFGGLYNDID DEFLVGTRKF IIRMRLSKLT YRDLVTNNII LMKRLKDIGH ISPEICRKYG ATGPVARGSG INFVVRKNEP YSVYNEFDFFD IPVYHEGDSY ARYMVRMDEM EQSLRIIEQA MDKLPGGPII PEKKPKIIKL PEGDYYSTVE AARGSGIRL VSDGGKNAYR LKLRSPTFNS MHLFDEVCQG MLIADALALM GSLDLVIPEI DR
Source	Baculovirus
Target Names	nuoD
Protein Names	Recommended name: NADH-quinone oxidoreductase subunit D EC= 1.6.99.5 Alternative name(s): NADH dehydrogenase I subunit D NDH-1 subunit D
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.