



Recombinant Rat Heterogeneous nuclear ribonucleoprotein H2 (Hnrnph2)

Product Code	CSB-EP731817RA-B
Abbreviation	Hnrnph2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AY09
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MMLSTEGREG FVVKVRGLPW SCSAEEVMRF FSDCKIQNGT SGVRFIYTRE GRPSGEAFVE LESEDEVKLA LKKDRETMGH RYVEVFKSNS VEMDWVLKHT GPNSPDTAND GFVRLRGLPF GCSKEEIVQF FSGLEIVPNG MTLPVDFQGR STGEAFVQFA SQEIAEKALK KHKERIGHRY IEIFKSSRAE VRTHYDPPRK LMTMQRPGPY DRPGAGRGYN SIGRGAGFER MRRGAYGGGY GGYDDYGGYN DGYGFGSDRF GRDLNYCFSG MSDHRYGDGG SSFQSTTGHC VHMRLPYRA TENDIYNFFS PLNPMRVHIE IGPDGRVTGE ADVEFATHED AVAAMAKDKA NMQHRYVELF LNSTAGTSGG AYDHSYVELF LNSTAGASGG AYGSQMMGGM GLSNQSSYGG PASQQLSAGY GGGYGGQSSM SGYDQVLQEN SSDYQSNLA
Source	E.coli
Target Names	Hnrnph2
Protein Names	Recommended name: Heterogeneous nuclear ribonucleoprotein H2 Short name= hnRNP H2 Alternative name(s): Heterogeneous nuclear ribonucleoprotein H' Short name= hnRNP H'
Expression Region	1-449
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.