



Recombinant Rat Sorting nexin-17 (Snx17)

Product Code	CSB-BP731834RA
Abbreviation	Snx17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AYS6
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MHFSIPETES RSGDSGGSTY VAYNIHVNGV LHCRVRYSQL LGLHEQLRKE YGANVLPAPF PKKLFSLTPA EVEQRREQLE KYMQAVRQDP LLGSSETFNS FLRRAQQETQ QVPTEEVSL E VLLSDGQKVL VNVLTSDQTE DVLEAVAAKL DLPDDLIGYF SLFLVREKED GAFSFVRKLQ EFELPYVSVT SLRSQEYKIV LRKSYWDSAY DDDVMENRVG LNLLYAQTVS DIEHWILVT KEQHRQLKSL QEKVSKKEFL RLAQTLRHYG YLRFDACVAD FPEKDCPVVV SAGNSELSLQ LRLPGQQLRE GSFVRTRMRC WRVTSSVPLP SGGTSTPSRG RGEVRLELAF EYLMKDRQLQ WVTITSPQAI MMSICLQSMV DELMVKKSSG SIRKMLRRRV GGNLRRSDSQ QAVKSPPLLE SPDASRESMV KLSSKLSAVS LRIGIGPSTD ASASAVHGNF AFEGIGDEDL
Source	Baculovirus
Target Names	Snx17
Protein Names	Recommended name: Sorting nexin-17
Expression Region	1-470
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein does not contain a coiled coil region, like some family members, but contains a B41 domain. This protein interacts with the cytoplasmic domain of P-selectin, and may function in the intracellular trafficking of P-selectin.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.