



Recombinant Debaryomyces hansenii H/ACA ribonucleoprotein complex subunit 3 (NOP10)

Product Code	CSB-BP731907DIS
Abbreviation	NOP10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BQP3
Product Type	Recombinant Protein
Immunogen Species	Debaryomyces hansenii (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (Torulaspora hansenii)
Purity	≥85% (SDS-PAGE)
Sequence	MHLMYTLGPD GKRVTYLNKT TSDGEITKSA HPARFSPDDK YSRQRVTLKK RFGLLPQTQ
Source	Baculovirus
Target Names	NOP10
Protein Names	Recommended name: H/ACA ribonucleoprotein complex subunit 3 Alternative name(s): Nucleolar protein 10 Nucleolar protein family A member 3 snoRNP protein NOP10
Expression Region	1-58
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.