



Recombinant *Yarrowia lipolytica* Pre-mRNA-splicing factor CWC21 (CWC21)

Product Code	CSB-EP731946YAP-B
Abbreviation	CWC21
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6C0M9
Product Type	Recombinant Protein
Immunogen Species	<i>Yarrowia lipolytica</i> (strain CLIB 122 / E 150) (Yeast) (<i>Candida lipolytica</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MSYNGIGLST PRGSATSGHI QTNISNRAFQ STNHRGEFAH HEMTDEQDKR LNKKLDRDYA ADRQADVELL EHERKRKVEV ACMELQDKLE EEGQKGEEI EEEVNTLRRQ LTLKMARELD VRTTIVSADR EALRKNYHKA AVIKQQEMDR MRDAFDTTKN LDDRGRDRPR NRQRQNREKS PVRDTRRSRS ARERSPSRRG RRPRRYSVTE DVRTSEVISE RQGSGSPCRR ADSRSRSRSP ERGSRRKSYR ERSRSRSRER SERRGRRDSV REREQDRSRS PLPY
Source	<i>E.coli</i>
Target Names	CWC21
Protein Names	Recommended name: Pre-mRNA-splicing factor CWC21
Expression Region	1-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.