



Recombinant Kluyveromyces lactis RuvB-like helicase 2 (RVB2)

Product Code	CSB-EP732067KBK
Abbreviation	RVB2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6CT29
Product Type	Recombinant Protein
Immunogen Species	Kluyveromyces lactis (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) (Candida sphaerica)
Purity	≥85% (SDS-PAGE)
Sequence	MSIQTQDSQE TLKSLSLIAA HSHISGLGLD ENLQPKPSSQ GMVGQLQARR AAGVILKMVQ NGSIAGRAVL VAGPPSTGKT ALAMGLSQSL GKDVPTAIA GSEIFSLELS KTEALTQAFR KSIGIKIEE TELIEGEVVE IQIDRSITGG HKQGKLTIKT TDMETIYELG NKMIEGLTKE KVLADGVDISI DKASGKITKL GRSFARSRDY DAMGADTKFV QCPEGELQKR KTVVHTVSLH EIDVINSRTQ GFLALFTGDT GEIRSEVRDQ INTKVAEWKE EGKAEIVPGV LFIDEVHMLD IECFSFINRA LEDEFAPIVI MATNRGISKT RGTTYKSPHG LPLDLLDRSI IITTKNYDAS EIKTILTIRS TEEVEELSPE ALDLLTNIGS ETSRLYSSNL ISVSQQAQK RKSNTVEVKD VERAYLLFLD SARSVKFVQE FESQYIDDKG QVNISIGKDE DAMDTTA
Source	E.coli
Target Names	RVB2
Protein Names	Recommended name: RuvB-like helicase 2 EC= 3.6.4.12
Expression Region	1-467
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.