



# Recombinant *Kluyveromyces lactis* RuvB-like helicase 2 (RVB2)

<b>Product Code</b>	CSB-YP732067KBK
<b>Abbreviation</b>	RVB2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6CT29
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Kluyveromyces lactis</i> (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) ( <i>Candida sphaerica</i> )
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSIQTQDSQE TLKSLSLIAA HSHISGLGLD ENLQPKPSSQ GMVGQLQARR AAGVILKMVQ NGSIAGRAVL VAGPPSTGKT ALAMGLSQSL GKDVPTAIA GSEIFSLELS KTEALTQAFR KSIGIKIKEE TELIEGEVVE IQIDRSITGG HKQGKLTIKT TDMETIYELG NKMIEGLTKE KVLADGVDISI DKASGKITKL GRSFARSRDY DAMGADTKFV QCPEGELQKR KTVVHTVSLH EIDVINSRTQ GFLALFTGDT GEIRSEVRDQ INTKVAEWKE EGKAEIVPGV LFIDEVHMLD IECFSFINRA LEDEFAPIVI MATNRGISKI RGTYNYKSPHG LPLDLLDRSI IITTKNYDAS EIKTILTIRS TEELEVELSPE ALDLLTNIGS ETSRLYSSNL ISVSQQIAQK RKSNTVEVKD VERAYLLFLD SARSVKFVQE FESQYIDDKG QVNISIGKDE DAMDTTA
<b>Source</b>	Yeast
<b>Target Names</b>	RVB2
<b>Protein Names</b>	Recommended name: RuvB-like helicase 2 EC= 3.6.4.12
<b>Expression Region</b>	1-467
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.