



Recombinant Kluyveromyces lactis Sterol 24-C-methyltransferase (ERG6)

Product Code	CSB-MP732096KKBK
Abbreviation	ERG6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6CYB3
Product Type	Recombinant Protein
Immunogen Species	Kluyveromyces lactis (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) (Candida sphaerica)
Purity	≥85% (SDS-PAGE)
Sequence	MSESNELRKR QAEFTKELHG NNVTKSGFSA LVSKNKSAQT EAVAKYLKHW DGKTDADAER RRLEDYNST HSYNVVTDF YEYGWGSSFH FSRFFKGESF SASVARHEHY LAYKAGIKEN DLILDVGCGV GGPARTISRF TGCNIIGLNN NDYQIQKANY YAKRDHLDK LSFVKGDFMK MEFDENTFDS VYAIEATCHA PTFEGVYGEI YKVLKPGGTF AVYEWVMTDK YDETNPENHRK IAYEIELGDG IPKMYSVDAA RDALSKVGFE ILEENDLADN DDEIPWYAPL TGEWKYVNSV ADLATFFRTS RLGRAFTTSM VTVFEKLGLA PKGSVSVTNA LEEAAVGLVA GGEKKLFTPM MLFVAKKPEQ K
Source	Mammalian cell
Target Names	ERG6
Protein Names	Recommended name: Sterol 24-C-methyltransferase EC= 2.1.1.41 Alternative name(s): Delta(24)-sterol C-methyltransferase
Expression Region	1-371
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.