



Recombinant *Erwinia carotovora* subsp. atroseptica Peptide methionine sulfoxide reductase MsrA (msrA)

Product Code	CSB-BP732121EAAB
Abbreviation	msrA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6D151
Product Type	Recombinant Protein
Immunogen Species	<i>Pectobacterium atrosepticum</i> (strain SCRI 1043 / ATCC BAA-672) (<i>Erwinia carotovora</i> subsp. atroseptica)
Purity	≥85% (SDS-PAGE)
Sequence	MNSIDKTQRI TQSDALPGRS TPMPVARLHV VHEHSMTHVP DHMSVAIFAM GCFWGAERLF WQQPGIYSTA AGYIGGYTPN PTYREVCSGQ TDHAEAVRVV FDPATISYSQ LLQLFWENHD PAQGMRRGGD IGSQYRSAIY TLTPEQEQA QASLQRFQQA MRENSDGRAI STEIEPAGPF YYAEEDHQQY LHKNPNGYCG LGGIGVCLPP PR
Source	Baculovirus
Target Names	msrA
Protein Names	Recommended name: Peptide methionine sulfoxide reductase MsrA Short name= Protein-methionine-S-oxide reductase EC= 1.8.4.11 Alternative name(s): Peptide-methionine (S)-S-oxide reductase Short name= Peptide Met(O) reductase
Expression Region	1-212
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.