



Recombinant *Erwinia carotovora* subsp. atroseptica Hemin import ATP-binding protein HmuV (hmuV)

Product Code	CSB-BP732154EAAB
Abbreviation	hmuV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6D645
Product Type	Recombinant Protein
Immunogen Species	<i>Pectobacterium atrosepticum</i> (strain SCRI 1043 / ATCC BAA-672) (<i>Erwinia carotovora</i> subsp. <i>atroseptica</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MTDSALDRTL VARHLRFQTN GRYLTDDVSL ELNCGEIVAI IGPNGAGKST LLRLLTGylt PDDGECllAG QPFSHWQPSA LAKTRAVMRQ HSGMAFAFSV QDVVAMGRSP HGRYPKNDDV VQQVMAQTGC LELATRDYRH LSGGEQQRVQ LARVLAQLWH PEPTGWLFL DEPTSALDLY HQQHLLRLLK QLTREQPLAV CCVLHDLNLA ALYADRILL HEGKLVAQGS PADVLQAETL AHWYRADLSV GSHPDYAIPQ VYLRQ
Source	Baculovirus
Target Names	hmuV
Protein Names	Recommended name: Hemin import ATP-binding protein HmuV EC= 3.6.3.-
Expression Region	1-265
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.