



Recombinant *Xenopus laevis* WD repeat-containing protein 37 (wdr37)

Product Code	CSB-MP732203XBE
Abbreviation	wdr37
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DDF0
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	≥85% (SDS-PAGE)
Sequence	MPTESGSWAA ARQTKQKRKS HSLSIKRTNS SEQDRLGLQR DMLEGQDSKL PPSVRNTLLE LFGQIEREFE NLYLENLELR REIDTLNDRL AVEGQAIDGA ELSKGQMKTK ASHSTSQLSQ KLKTTYKAST SKIVSSFKTT TSRAICQLVK DYVGHRDGLW DVSVRTQPV VLGTASADHT ALLWSIETGK CLIKYVGHAG SVNSIKFHPT EQIALTASGD QTAHIWRYMV QLPTPQPMAD TSQISGEEEV DFSDKDENDG DGDASSDCPT VRVPLTALKS HQGVVIAADW LVGGKQAVTA SWDRTANLYD VETSELVHSL TGHDQELTHC CTHPTQRLVV TSSRDITFRL WDFRDPSIHS VNVFQGHTDT VTSAVFTVGD NVVSGSDDRT VKVWDLKNMR SPIATIRTDS AVNRISVSVG QRIIALPHDN RQVRLFDISG VRLARLPRSN RQGHRRMVCC CAWSEDHPTC NLFSCGFDRQ AIGWNINIPA LLQEK
Source	Mammalian cell
Target Names	wdr37
Protein Names	Recommended name: WD repeat-containing protein 37
Expression Region	1-495
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.