



Recombinant Danio rerio Coiled-coil domain-containing protein 75 (ccdc75)

Product Code	CSB-EP732231DIL
Abbreviation	ccdc75
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DGZ0
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MADDEDDYMS DAFLNQISDV RPGVPMVKRV KEAIKKEALH KEKDIQNRQK SFKEQEKESR EAALQSSIGS QNKGFALLQK MGYKAGQGLG KEGAGRVEPV PLNIKTDRGG IGMEEVKKRK ADEELQNYRR KVHMKQHLEK KSIEDFRVRK RTEREERQTQ GDLRKSQRAC EQLDNQKGIT VPRDCWYWPE VEKNEESELL EEEIKTEESD EEEELTPLDK LQFLTSYLRG VHYYCIWCGT AYNDEEDLGT NCPGDTAADH DD
Source	E.coli
Target Names	gpatch11
Protein Names	Recommended name: Coiled-coil domain-containing protein 75
Expression Region	1-262
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.