



Recombinant Danio rerio Coiled-coil domain-containing protein 85C-B (ccdc85cb)

Product Code	CSB-MP732235DIL
Abbreviation	ccdc85cb
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DHL7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MAKHCPNDDL SKIPDEELVR WSKEDLIRRL RRV DGEKMSL MLEHSNMMKD VNRSLQVHLH EIRNLKEINQ KLQDDNQELR ELCCFLDDDR QKGKKSREW QRF GKFTASS VWKEVSTYQQ KLKELEINQE NVLRENAELK DIILMLDEDR NGAGSRSSID SQSSLSNLNG GSGTVRDVGD GSSTSSGGSA GSPDHHHNI HKTVEAKIGT VRRSMDDLSA PHHHR SIPAG LNDASSNYIR QLETKVRILE DNNKLLSQP CSRYLSKLS MKYNPGDLRA LRKGMTLYHS ESQ LSSLPQR QEALLNGTGR LQTSESSPST GFISSAQKPE AVVHAMKVLE VHDNLEKQLP EESEEDLSEK EKAIVREMCN VVWRKLG DAT GTKPSLRQQL SGNQFKAPL
Source	Mammalian cell
Target Names	ccdc85cb
Protein Names	Recommended name: Coiled-coil domain-containing protein 85C-B
Expression Region	1-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.