



Recombinant Human TBC1 domain family member 3G (TBC1D3G)

Product Code	CSB-MP732237HU
Abbreviation	TBC1D3G
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DHY5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDVVEVAGSW WAQEREDIIM KYEKGHRAGL PEDKGPKPFR SYNNNVVDHLG IVHETELPPL TAREAKQIRR EISRKSKWVD MLGDWEKYKS SRKLIDRAYK GMPMNIRGPM WSVLLNIEEM KLKNPGRYQI MKEKGKRSSE HIQRIDRDIS GTLRKHMFFR DRYGTKQREL LHILLAYEEY NPEVGYCRDL SHIAALFLLY LPEEDAFWAL VQLLASERHS LQGFHSPNGG TVQGLQDQQE HVVATSQPKT MGHQDKKDLG GQCSPLGCLI RILIDGISLG LTLRLWDVYL VEGEQALMPI TRIAFKVQQK RLTKTSRCGP WARFCNRFVD TWARDLTVL KHLRASMKKL TRKQGDLPPL AKPEQGSSAS RPVPASRGGK TLCKGDRQAP PGPPARFPRP IWSASPPRAP RSSTPCPGGA VREDTYPVGT QGVPSPALAQ GGPQGSWRFL QWNSMPRLPT DLDVEGPWFR HYDFRQSCWV RAISQEDQLA PCWQAEHPAE RVRSAFAAPS TDSDQGTPFR ARDEQQCAPT SGPCLCGLHL ESSQFPPGF
Source	Mammalian cell
Target Names	TBC1D3G
Protein Names	Recommended name: TBC1 domain family member 3G
Expression Region	1-549
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.