



Recombinant *Xenopus laevis* Protein CWC15 homolog B (cwc15-b)

Product Code	CSB-BP732252XBE
Abbreviation	cwc15-b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DKE6
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	≥85% (SDS-PAGE)
Sequence	MTTAARPTFE PARGGRGKGE SDLSQLSKQY SSRDLPSHTK IKYRQTTQDA PEEVRSRDFR RELEERERVA VRDKNRDRPT REHTSSVSKK PRLDQIPAAN LDADDPLTDE DADEDSDEDS DDDTAALLAE LEKIKKERAE EQVRKELEQK AEEERIRMEN ILSGNPLLNL TGPAAQSQAS FKVKRRWDDD VVFKNCAKGV DEMKKRQDKR FVNDTLRSEF HKKFMEKYVK
Source	Baculovirus
Target Names	cwc15-b
Protein Names	Recommended name: Protein CWC15 homolog B
Expression Region	1-230
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.