



Recombinant Gecko japonicus Protein TSSC1 (TSSC1)

Product Code	CSB-EP732272GEG-B
Abbreviation	TSSC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DUZ9
Product Type	Recombinant Protein
Immunogen Species	Gekko japonicus (Schlegel's Japanese gecko)
Purity	≥85% (SDS-PAGE)
Sequence	MEDDAPVIYG LEFQARALTP QTAETDAIRF LVGTQSLKYD NQIHIDFDD ENNIINKNVL LHQVGEIWHI SASPADKGVL ATCYSKTSDS KVMTC AAVWR MPKELESGSH ESPDDSSSNA QTLELLCHLD NTAHGNMAGV LWPEMGDGKK LISLADNHLL LWDLQESSK AVLSNSAALE GKGQLKFTSG RWSPHHNCTQ VATANDTAIR GWDIRTMSQI YCIESAHGQL VRDLDFNPNK QYYLASCDD CKVKFDWTRN VSDPVKTL EE HSHWVWNVRY NHSHDQLVLT ASSDSRVILS NMVSISEPF GHLVDDEELS DQEDQH QEDK IKEPLQDSII ATYEEHEDSV YAVEWSSADP WLFASLSYDG RLVINRV PRA LKYHILL
Source	E.coli
Target Names	eipr1
Protein Names	Recommended name: Protein TSSC1
Expression Region	1-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.