



Recombinant Human Putative UPF0607 protein LOC392364

Product Code	CSB-EP732509HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PR19
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAPAPVPALG TLGCYRFLFN PRQHLGPSFP ARRYGAPRRL CFLPQNTGTP LRVLPSVFWS PPSRKKPVL ARNSRMFGHL SPVRIPHLRG KFNLRPLSLD EQVIPARLPK MEVRAEEPKE ATEVKDQVET QEVEDNKRGP CSNGEAASTS RPLETQGNPT SPRYNPRPLE GNVQLKSLTE NNQTDKAQVH AVSFYSKGHG VASSHSPAGG FPRGRTPPAR QHG
Source	E.coli
Protein Names	Recommended name: Putative UPF0607 protein LOC392364
Expression Region	1-223
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.