



Recombinant Danio rerio Eukaryotic translation initiation factor 3 subunit H-B (eif3hb)

Product Code	CSB-EP732510DIL-B
Abbreviation	eif3hb
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PR67
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MASRKETPNG KTLDTVPVKQI QIEGLVVMKI IKHYQEEGQG SEVVQGVLLG LVVEDQLEIT NCFPPQHT E DDADFDEVQY QMEMMRSLRH VNIDHLHVGW YQSTLYGSFV SRALLDSQFS YQHAIEESVV LIYDPIKTAQ GSLCLKAYRL TPKLMEICKE KDFSSEGLKK ASVGYEHMFE EVPIVIKNSH LINVLLWELE EKSTAADKHE LLSLSSSSHL EKSLQMLMDR VDDMSQDIVK YNNYSRSLSK QQQQKHQYVQ RRQQENAQRQ SRGEPPLPEE DLTKMFKPPQ PPPRMDTLI ASQINTYCQT IKEFTSQNLG KLFMAEALQG PSS
Source	E.coli
Target Names	eif3hb
Protein Names	Recommended name: Eukaryotic translation initiation factor 3 subunit H-B Short name= eIF3h-B Alternative name(s): Eukaryotic translation initiation factor 3 subunit 3-B eIF-3-gamma-B eIF3 p40 subunit B
Expression Region	1-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.