



# Recombinant Pongo abelii Dematin (EPB49)

<b>Product Code</b>	CSB-MP732593PYX
<b>Abbreviation</b>	EPB49
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5R4B6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MERLQKQPLT SPGSVSPSRD SSVPGSPSSI VAKMDNQVLG YKDLAaipkd KAILDIERPD LMIYEPHFTY SLLEHVELPR SRERSLSPKS TSPPPSPEVW ADSRSPGIIS QASAPRTTGT PRTSLPHFHH PETS RPDSNI YKKPPIYKQR ESVGGSPQTK HLIEDLIIES SKFPAAQPPD PNQPAKIETD YWPCPPSLAV VETEWKRKA SRRGAE EEEEE EEDDDSGEEM KALRERQREE LSKVTSNLGK MILKEEMEKSLPIRRKTRSL PDRTPFHTSL HQGTSKSSSL PAYGRTTSLR LQSTEFSPSG SETGSPGLQN GEGQRGRMDR GNSLPCVLEQ KIIPPYEMLVV TNKGRTKLPP GVDRMLERH LSAEDFSRVF AMSPEEFGKL ALWKRNELKK KASLF
<b>Source</b>	Mammalian cell
<b>Target Names</b>	DMTN
<b>Protein Names</b>	Recommended name: Dematin Alternative name(s): Erythrocyte membrane protein band 4.9
<b>Expression Region</b>	1-405
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.