



Recombinant Pongo abelii Vasculin (GPBP1)

Product Code	CSB-MP732595PYX
Abbreviation	GPBP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R4E9
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MAQHDFAPAW LNFPTPPSST KSSLNFEKHS ENFAWTENRY DVNRRRRHNSS DGFDSAIGRP NGGNFGRKEK NGWRTHGRNG TENINHRGGY HGGSSRSRSS IFHSGKSQGL HENNIPDNET GRKEDKRERK QFEAEDFPSL NPEYEREPNH NKSLAAGVWE YPPNPKSRAP RMLVIKKGNT KDLQLSGFPV VGNLPSQPVK NGTGPSVYKG LVPKPAAPPT KPTQWKSQTK ENKVGTSFPH ESTFGVGNFN AFKSTAKNFS PSTNSVKECN RSNSSSPVDK LNQQPRLTKL TRMRTDKKSE FLKALKRDRV EEEHEDESRA GSEKDDDSFN LHNSNSTHQE RDINRNF DEN EIPQENG NAS VISQQIIRSS TFPQTDVLSS SLEAEHRLK EMGWQEDSEN DETCAPITED EMREFQVISE QLQKNGLRKN GILKNGLICD FKFGPWKNST FKPTTENDDT ETSSSDTSDD DDV
Source	Mammalian cell
Target Names	GPBP1
Protein Names	Recommended name: Vasculin Alternative name(s): GC-rich promoter-binding protein 1
Expression Region	1-473
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.