



# Recombinant Pongo abelii Mitochondrial-processing peptidase subunit alpha (PMPCA)

<b>Product Code</b>	CSB-MP732607PYX
<b>Abbreviation</b>	PMPCA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5R513
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SSGGAYP NIPLSSPLPG VPKPVFATVD GQEKFKTKVT TLDNGLRVAS QNKFGQFCTV GILINSGSRY EAKYLSGIAH FLEKLAFSST ARFDSKDEIL LTLEKHGGIC DCQTSRDTTM YAVSADSKGL DTVVGLLADV VLQPRLTDEE VEMTRMTVQF ELEDNLNLRPD PEPLLTEMIH EAAYRENTVG LHRFCPTENI AKINREVLHS YLRNYYTPDR MVLAVGVGEH EHLVDCARKY LLGIQPAWGS AEAVDIDRSV AQYTGGIKR ERDMSNVSLG PTPIPELTHI MVGLESCSFL EEDFIPFAVL NMMMGSGGSF SAGGPGKGMF SRLYLNLVLR HHWMYNATSY HHSYEDTGLL CIHASADPRQ VREMVEITK EFILMSGTVD AVELERAKTQ LTSMLMMNLE SRPVIFEDVG RQVLATRSRK LPHELCTLIR NVKPEDVKRV ASKMLRGKPA VAALGDLTDL PTYEHIQTAL SSKDGRLPRT YRLFR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PMPCA
<b>Protein Names</b>	Recommended name: Mitochondrial-processing peptidase subunit alpha EC=3.4.24.64 Alternative name(s): Alpha-MPP
<b>Expression Region</b>	34-525
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.