



# Recombinant Pongo abelii Methyltransferase-like protein 14 (METTL14)

<b>Product Code</b>	CSB-MP732620PYX
<b>Abbreviation</b>	METTL14
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5R5N4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDSRLQEIRE RQKLRRQLLA QQLGAESADS IGAVLNSKDE QREIAETRET CRASYDTSAP NAKRKYLDEG ETDEDKMEFY KDELEMQQDE ENLPYEEEEIY KDSSTFLKGT QSLNPHNDYC QHFVDTGHRP QNFIRDVGLA DRFEEYPKLR ELIRLKDELI AKSNTPPMYL QADIEAFDIR ELTPKFDVIL LEPPLEYYR ETGITANEKC WTWDDIMKLE IDEIAAPRSF IFLWCGSGEG LDLGRVCLRK WGYRRCEDIC WIKTNKNNPG KTKTLDPAKAV FQRTKEHCLM GIKGTVKRST DGDFIHANVD IDLIITEEPE IGNIKPVFI FHIIEHFCLG RRRLHLFGRD STIRPGWLV GPTLTNSNYN AETYASYFSA PNSYLTGCTE EIERLRPKSP PPKSKSDRGG GAPRGGGRRG TSAGRGRERN RSNFRGERGG FRGGRGGAHR GGFPPR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	METTL14
<b>Protein Names</b>	Recommended name: Methyltransferase-like protein 14 EC= 2.1.1.-
<b>Expression Region</b>	1-456
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.