



# Recombinant Pongo abelii Cyclin-dependent kinase-like 2 (CDKL2)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-BP732647PYX                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Abbreviation</b>      | CDKL2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                              |
| <b>Uniprot No.</b>       | Q5R754                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Immunogen Species</b> | Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Sequence</b>          | MEKYENLGLV GEGSYGMVMK CRNKDTGRIV AIKKFLESDD DKMVKKIAMR<br>EIKLLKQLRH ENLVNLLVC KKKRWYLVF EFVDHTILDD LELFPNGLDY<br>QVVQKYLFQI INGIGFCHSH NIIHRDIKPE NILVSKSGVV KLCDFGFART<br>LAAPGEVYTD YVATRWRAP ELLVGDKYK KAVDVWAIGC LVTEMFMGEP<br>LFPGDSIDQ LYHIMMCLGN LIPRHQELFY KNPVFAGVRL PEIKEREPL<br>RRYPKLSEVV IDLAKKCLHI DPKRPFCAE LLHHDFQMD GFAERFSQEL<br>QLKVQKDARN VLSKKSQNR KKEKEKDDSL GEERKTLVVQ DTNADPKIKD<br>CKLFIKIGSK IDGEKAEKGS RASNASCLHD SRTSHNKIVP STSLKDCSNV<br>SVDHTRNPSV AIPPLTHNLS AVAPGINSGM GTATPIQGH RVDEKTKKCS<br>IPFVKPNRHS PSGIYNINVT TLVSGPPLSD DSGADLPQME HQH |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Target Names</b>      | CDKL2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Protein Names</b>     | Recommended name: Cyclin-dependent kinase-like 2 EC= 2.7.11.22                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Expression Region</b> | 1-493                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                     |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life                                                                                                                                                                                                                                                                                                                                               |



of lyophilized form is 12 months at -20°C/-80°C.