



Recombinant Pongo abelii KAT8 regulatory NSL complex subunit 2 (KANSL2)

Product Code	CSB-MP732666PYX
Abbreviation	KANSL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R802
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MNRIRIHVLP TNRGRITPVP RSQEPLSCAF THRPCSQPRL EGQEFCIKHI LEDKDAPFKQ CSYISTKNGK RCPNAAPKPE KKDGVSFCAE HVRRNALALH AQMKKTNPGP MGETLLCQLS SYAKTELGSQ TPRESSRSEAS RILDEDSWSD GEQEPITVDQ TWRGDPDSEA DSIDSDQEDP LKHAGVYTAE EVALIMREKL IRLQSLYIDQ FKRLQHLLKE KKRRYLHNRK VEHEALGSSL LTGPEGLLAK ERENLKRLKC LRRYRQRYGV EALLHRQLKE RRMLATDGAA QQAHTTRSSQ RCLAFVDDVR CSNQSLPMTR HCLTHICQDT NQVLFKCCQG SEEVPCNKPV PVSLSEDPCC PLHFQLPPQM YKPEQVLSVP DDLEAGPMDL YLSAAELQPT ESLPLEFSDD LDVVGDMQC PPSPLLFDPs LTLEDHLVKE IAEDPVDILG QMQMAGDGCR SQGSRNSEKG SAPLSQSGLA TANGKPEPTS IS
Source	Mammalian cell
Target Names	KANSL2
Protein Names	Recommended name: KAT8 regulatory NSL complex subunit 2 Alternative name(s): NSL complex protein NSL2 Non-specific lethal 2 homolog
Expression Region	1-492
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.