



Recombinant Pongo abelii WD repeat-containing protein 89 (WDR89)

Product Code	CSB-EP732746PYX-B
Abbreviation	WDR89
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RBZ3
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MEKIEEQFAN LHIVKCSLGT KEPTYLLGID TSKTVQAGKE NLVAVLCSNG SIRIYDKERL SVLREFSGYP GLLNGVRFAN SCDSVYSACT DGTVKCWDAR VAREKPVQLF KGYP SNIFIS FDINCNDHII CAGTEKVDDD ALLVFW DARM NSQDLSTTKD PLGAYSETHS DDVTQVRFHP SNP NMVVS GS SDGLVNVFDI NIDNEEDALV TTCNSISSVS CIGWSGKGYK QIYCMTHDEG FYWWDLNHLD TDEPVTRLNI QDVREVVNMK EDALDYLIGG LYHEKTD TLH VIGGTNKGRI HLMNCSV SGL THV TSLQGGH AATVRSFCWN VQDDSLLTGG EDAQLLLWKP GAIEKTFTKK ESMKIASSVH QRVRVHSNDS YKRRKKQ
Source	E.coli
Target Names	WDR89
Protein Names	Recommended name: WD repeat-containing protein 89
Expression Region	1-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.