



# Recombinant Pongo abelii WD repeat-containing protein 89 (WDR89)

<b>Product Code</b>	CSB-MP732746PYX
<b>Abbreviation</b>	WDR89
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5RBZ3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEKIEEQFAN LHIVKCSLGT KEPTYLLGID TSKTVQAGKE NLVAVLCSNG SIRIYDKERL SVLREFSGYP GLLNGVRFAN SCDSVYSACT DGTVKCWDAR VAREKPVQLF KGYP SNIFIS FDINCNDHII CAGTEKVDDD ALLVFW DARM NSQDLSTTKD PLGAYSETHS DDVTQVR FHP SNP NMVVS GS SDGLVNVFDI NIDNEEDALV TTCNSISSVS CIGWSGKGYK QIYCMTHDEG FYWWDLNHLD TDEPVTRLNI QDVREVVNMK EDALDY LIGG LYHEKTD TLH VIGGTNKGRI HLMNCSV SGL THV TSLQGGH AATVRSFCWN VQDDSLLTGG EDAQLLLWKP GAIEKTFTKK ESMKIASSVH QRVRVHSNDS YKRRKKQ
<b>Source</b>	Mammalian cell
<b>Target Names</b>	WDR89
<b>Protein Names</b>	Recommended name: WD repeat-containing protein 89
<b>Expression Region</b>	1-387
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.