



Recombinant Thermus thermophilus Putative regulator of ribonuclease activity (TTHA1322)

Product Code	CSB-MP732873TNT
Abbreviation	TTHA1322
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SIP7
Product Type	Recombinant Protein
Immunogen Species	Thermus thermophilus (strain HB8 / ATCC 27634 / DSM 579)
Purity	≥85% (SDS-PAGE)
Sequence	MEARTTDLSD LYPEGEALPM VFKSFGGRAR FAGRVRTLRLV FEDNALVRKV LEEEGAGQVL FVDGGGSLRT ALLGGNLARL AWEKGWAGVV VHGAVRDTEE LREVPIGLLA LAATPKKSAK EGKGEVDVPL KVLGVEVLPG SFLLADEDGL LLLPEPPSGV RSGG
Source	Mammalian cell
Target Names	TTHA1322
Protein Names	Recommended name: Putative regulator of ribonuclease activity
Expression Region	1-164
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.