



Recombinant *Thermus thermophilus* Ribosomal RNA small subunit methyltransferase H (rsmH)

Product Code	CSB-BP732879TNT
Abbreviation	rsmH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SJD8
Product Type	Recombinant Protein
Immunogen Species	<i>Thermus thermophilus</i> (strain HB8 / ATCC 27634 / DSM 579)
Purity	≥85% (SDS-PAGE)
Sequence	MRPMTHVPVL YQEALDLLAV RPPGGVYVDAT LGGAGHARGI LERGGRVIGL DQDPEAVARA KGLHLPGLTV VQGNFRHLKR HLAALGVERV DGILADLGVS SFHLDDPSRG FSYQKEGPLD MRMGLEGPTA KEVVNRLPLE ALARLLRELG EEPQAYRIAR AIVAAREKAP IETTTQLAEI VRKAVGFERRA GHPARKTFQA LRIYVNDELN ALKEFLEQAA EVLAPGGRLV VIAFHSLEDR VVKRFLRESG LKVLTKKPLV PSEKEAAQNP RARSACLRAA EKEAP
Source	Baculovirus
Target Names	rsmH
Protein Names	Recommended name: Ribosomal RNA small subunit methyltransferase H EC=2.1.1.199 Alternative name(s): 16S rRNA m(4)C1402 methyltransferase rRNA (cytosine-N(4)-)-methyltransferase RsmH
Expression Region	1-285
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.