



Recombinant Mouse Protein FAM71B (Fam71b)

Product Code	CSB-MP732906MO
Abbreviation	Fam71b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5STT6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MPGMKRTMSS ECLLPYYTAH SYRSMGVFNT SMGNLQRQLY KGGEYDIFKY APMFESDFIQ ISKRGEVIDV HNRVRMVTVC IASTSPVLPL PDVMLLARPA KVCEEHARRA RFIKKGKCKP SKTLELTRLL PLKFKVISVH DHEKQQLRLK LATGRTFYLQ LCPSSDARED LFCYWEKLVY LLRPPMTNCI SNSTLPTGET SVDTKSTLVS EIRGEGDQNS RPQSSPTVSE ATSAAFAGGE RTQPAAAAAV TPVSAKARAA GTAGAAAGTA GAAAGTAGPA AGPAAGTAGP AAGTAGAAAG TAGAAAGTAG ATAGTAGATA GMAGATAGTA AETAGTAAET AGAARAAGGP MAAAVAAPSA GMTKAETATS PTSGVISLAA TTTKPPGSGQ VAAAMIGSAA KDQVGGESSK AMALANITL ENVDVALAGA ANSISESPPA GGDASGSPDT GLNVAFAGSI KTKSPAEDKP EAPLVSTLQS EGYMCERDGS QKVSQTSSEA KEKRERREKD RTSSRKSSH RRTGMSRHSS KDKSRKTSSY RSVSGKTRED KGKGHGRLRG KRHSSSHKSE SRTGHKTRKN RSPAGLGSVS KRATKITSFF RSFLVRPTPK AGDTSCDRGG VDIVTKLVEK KQDIEAVMEK SKDSEFKDTV ISETMEKIIL ETKSI
Source	Mammalian cell
Target Names	Fam71b
Protein Names	Recommended name: Protein FAM71B
Expression Region	1-665
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.