



# Recombinant Mouse Protein FAM71B (Fam71b)

<b>Product Code</b>	CSB-YP732906MO
<b>Abbreviation</b>	Fam71b
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5STT6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPGMKRTMSS ECLLPYYTAH SYRSMGVFNT SMGNLQRQLY KGGEYDIFKY APMFESDFIQ ISKRGEVIDV HNRVRMVTVC IASTSPVLPL PDVMLLARPA KVCEEHARRA RFIKKGKCKP SKTLELTRLL PLKFKVISVH DHEKQQLRLK LATGRTFYLQ LCPSSDARED LFCYWEKLVY LLRPPMTNCI SNSTLPTGET SVDTKSTLVS EIRGEGDQNS RPQSSPTVSE ATSAAFAGGE RTQPAAAAAV TPVSAKARAA GTAGAAAGTA GAAAGTAGPA AGPAAGTAGP AAGTAGAAAG TAGAAAGTAG ATAGTAGATA GMAGATAGTA AETAGTAAET AGAARAAGGP MAAAVAAPSA GMTKAETATS PTSGVISLAA TTTKPPGSGQ VAAAMIGSAA KDQVGGESSK AMALANITL ENVDVALAGA ANSISESPPA GGDASGSPDT GLNVAFAGSI KTKSPAEDKP EAPLVSTLQS EGYMCERDGS QKVSQTSSEA KEKRERREKD RTSSRKSSH RRTGMSRHSS KDKSRKTSSY RSVSGKTRED KGKGHGRLRG KRHSSSHKSE SRTGHKTRKN RSPAGLGSVS KRATKITSFF RSFLVRPTPK AGDTSCDRGG VDIVTKLVEK KQDIEAVMEK SKDSEFKDTV ISETMEKIIL ETKSI
<b>Source</b>	Yeast
<b>Target Names</b>	Fam71b
<b>Protein Names</b>	Recommended name: Protein FAM71B
<b>Expression Region</b>	1-665
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.