



Recombinant Mouse DnaJ homolog subfamily C member 11 (Dnajc11)

Product Code	CSB-EP732970MO-B
Abbreviation	Dnajc11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5U458
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	ATALSEEEL DNEDYYSLLN VRREASSEEL KAAAYRRLCML YHPDKHRDPE LKSQAERLFN LVHQAYEVLS DPQTRAIYDI YGKRGLEMEG WEVVERKRTP AEIREEFERL QREREERRLQ QRTNPKGTIS VGVDATDLFD RYDEEYEDVS GSGFPQIEIN KMHISQSIEA PLTATDTAIL SGSLSTQNGN GGGSVNFALR RVTSAGWGE LEFGAGDLQG PLFGLKLFNRN LTPRCFVTTN CALQFSSRGI RPGLTTVLAR NLDKNTVGYL QWRWGIQSAM NTSIVRDTKT CHFTVALQLG IPHSFALISY QHKFQDDDQT RVKGSCLKAGF FGTIVEYGAE RKISRHSVLG AAVSIGVPQG VSLKVKNRA SQTYFFPIHL TDQLLPSAVF YATVGPLVVY LAVHRLIIRP YLRAQKEKEL EKQRENTASD ILQKKQEAEA AVRLMQESVR RIIEAEESRM GLIIVNAWYG KFVNDKSRKN EKVKVIDVTV PLQCLVKDSK LILTEASKAG LPGFYDPCVG EEKSLRVLYQ FRGVLHQVMV PDSEALRIPK QSHRIDTDG
Source	E.coli
Target Names	Dnajc11
Protein Names	Recommended name: DnaJ homolog subfamily C member 11
Expression Region	2-559
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.