



Recombinant Human Probable E3 ubiquitin-protein ligase HECTD2 (HECTD2)

Product Code	CSB-EP732977HU-B
Abbreviation	HECTD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5U5R9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSEAVRVPSP ATPLVVAAPA PEERKKGKES REKLPIVSA GAGATAGLDR GAKGQISTFS SFISAVSPKK EAAENRSSPA HLVPFNIKNV REPPPICLDV RQKQRTSMDA SSSEMKAAPVL PEPILPIQPK TVKDFQEDVE KVKSSGDWKA VHDFYLTTFD SFPELNAAFV KDATASFNTI EDSGINAKFV NAVYDTLLNT PQDVQKTVLK GIINSLLEW KGPRTKDDL R AYFILLQNPQ FNNTSTYVIY AHLRQIATL VEADHHFLVH WFKKLSQKRF KQLVERLLQF ISLRLFPAKP EEFPITKCS WWIPSAKVL ALLNTANNLV HPPLIPYTFD YNSTLDHIDL MEEYHTWQNF GNSHRFSFCQ YPFVISVAAK KIIQRDSEQ QMINIARQSL VDKVSRRQRP DMNILFLNMK VRRTHLVSDS LDELTRKRAD LKKKLKVTFF GEAGLDMGGL TKEWFLLIR QIFHPDYGMF TYHKDSHCHW FSSFKCDNYS EFRLVGILMG LAVYNSITLD IRFPCCYKK LLSPPIIPSD QNIPVGICNV TVDDLQIMP ELAHLSELL SHEGNVEEDF YSTFQVFQEE FGIKSYNLK PGGDKISVTN QNRKEYVQLY TDFLLNKSIV KQFAAFYYGF HSVCASNALM LLRPEEVEIL VCGSPDLDMH ALQRSTQYDG YAKTDLTIKY FWDVVLGFPL DLQKLLHFT TGSDRVPVGG MADLNFKISK NETSTNCLPV AHTCFNQLCL PPYKSKKDLK QKLIIGISNS EGFGLE
Source	E.coli
Target Names	HECTD2
Protein Names	Recommended name: Probable E3 ubiquitin-protein ligase HECTD2 EC= 6.3.2.- Alternative name(s): HECT domain-containing protein 2
Expression Region	1-776
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.