



# Recombinant Human Probable E3 ubiquitin-protein ligase HECTD2 (HECTD2)

<b>Product Code</b>	CSB-MP732977HU
<b>Abbreviation</b>	HECTD2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5U5R9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSEAVRVPSP ATPLVVAAPA PEERKKGKES REKLPIVSA GAGATAGLDR GAKGQISTFS SFISAVSPKK EAAENRSSPA HLVFPNIKNV REPPPICLDV RQKQRTSMDA SSSEMKAAPVL PEPILPIQPK TVKDFQEDVE KVKSSGDWKA VHDFYLTTFD SFPELNAAFK KDATASFNTI EDSGINAKFV NAVYDTLLNT PQDVQKTVLK GIINSLLEW KGPRTKDDL R AYFILLQNPQ FNNTSTYVIY AHLRQIATL VEADHHFLVH WFKKLSQKRF KQLVERLLQF ISLRLFPAKP EEFPITKCS WWIPSAAKVL ALLNTANNLV HPPLIPYTFD YNSTLDHIDL MEEYHTWQNF GNSHRFSFCQ YPFVISVAAK KIIQRDSEQ QMINIARQSL VDKVSRRQRP DMNILFLNMK VRRTHLVSDS LDELTRKRAD LKKKLKVTFFV GEAGLDMGGL TKEWFLLIR QIFHPDYGMF TYHKDSHCHW FSSFKCDNYS EFRLVGILMG LAVYNSITLD IRFPCCYKK LLSPPIIPSD QNIPVGICNV TVDDLQIMP ELAHLSELL SHEGNVEEDF YSTFQVFQEE FGIKSYNLK PGGDKISVTN QNRKEYVQLY TDFLLNKSIIY KQFAAFYYGF HSVCASNALM LLRPEEVEIL VCGSPDLDMH ALQRSTQYDG YAKTDLTIKY FWDVVLGFPL DLQKLLHFT TGSDRVPVGG MADLNFKISK NETSTNCLPV AHTCFNQLCL PPYKSKKDLK QKLIIGISNS EGFGLE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	HECTD2
<b>Protein Names</b>	Recommended name: Probable E3 ubiquitin-protein ligase HECTD2 EC= 6.3.2.- Alternative name(s): HECT domain-containing protein 2
<b>Expression Region</b>	1-776
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.