



Recombinant Bombyx mori 60S ribosomal protein L17 (RpL17)

Product Code	CSB-MP732983BTT
Abbreviation	RpL17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UAS2
Product Type	Recombinant Protein
Immunogen Species	BOMMO(Silk moth)
Purity	≥85% (SDS-PAGE)
Sequence	MGARYSREPDN PAKSCKARGS NLRVHFKNTY ETAMAIRKMP LRRAVRYLKN VIEKKECIPF RRFNGGVGRC AQAKQFGTTQ GRWPKKSAEF LLQLLRNAES NADNKTLDVD RLVIDHIQVN RAPCLRRRTY RAHGRINPYM SSPCHIEVCL SEREDAVARV APTDDAPAKK KLSKKKLARQ KEKMMRE
Source	Mammalian cell
Target Names	RpL17
Protein Names	Recommended name: 60S ribosomal protein L17
Expression Region	1-187
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.