



Recombinant Acanthamoeba polyphaga mimivirus Uncharacterized protein R367 (MIMI_R367)

Product Code	CSB-BP733085ADAZ
Abbreviation	MIMI_R367
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UQV7
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MNCQNKILNS VTVVYGYSSG DFSWKNNDDR KIFRTSGLIA SIEGKKYVVT TREQMISCKT LIMYHLFKKN GVLVKSVMHI LFQSIENLI ILGSVGYNEL ILDKDNAVST TLNQKEIEII LNKFKGTNLD SHQTIPTNRS LYHINKIDLG NTFDNIKYKY HTYDIKFNKQ CIFDKSYLPK SLKYKFCQK NINSVGTGTLGT VVYSAKNKII GITTHVGTNK IYVTPTTVII SMMRNFIDYF NDQQSYKGQS CFPFKLIKIDE SNNIITMOK NKSNNKVPVK NGDMLRTLNG YEIIVNDNEA CLCVDSKNII PLDIYFRMKY QKDASYELTV SRGSKIVKFN LFLVSNTSEL RLTDQTYYP KDTIPHINLF GLIVVKLSHE LIDYLVKNGI IVRNYVIENI MNGAELNDDI LIILDSTNEL NCTTFNFPIL TKQTKPIDCP VLTSLNHQKI SNLENLSNVK LFTKGENILT LKLNQDSEID LLV
Source	Baculovirus
Target Names	MIMI_R367
Protein Names	Recommended name: Uncharacterized protein R367
Expression Region	1-483
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.