



Recombinant Acanthamoeba polyphaga mimivirus Uncharacterized protein R883 (MIMI_R883)

Product Code	CSB-YP733087ADAZ
Abbreviation	MIMI_R883
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UQX5
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MTTTNNVLLV DDQHRNIGVT DDLHTIESLK SELTLDYDLP TLCGKNLDSP NILLNSIEKF KETFNKHYPF LEKINMKNLL IAGGSVSNIV RNKFQYGSDI DFFIYGLNQQ EASSRVRQWL IDILVKKPDN SSEKTTTKKI DMTKYKIIIR NNNCITILLD HGDLKLQLIF RLYQSISEIL HGFDLGSSAV GYDGENVYFT TLGKYCHEYS CNVIDTTRS TTYEYRLNKY FDRGFNIVVP KLDLSKLTFT NLKYGEVEIC ELPYFIFSYQ NIIGNKIIVK KFYDKYNIKS DYGLEPINST NLYYQSLKIN IGNLINDVDY YYYVSSHIDK QNTDILTKAP RLTKGDIITF YDGVRTKLNG KNIDVCLIRK YITIDTIENI VATMFHKDTN VKEYFDSIIE KQKELALNKL DILLQKNHNI VWITENPGKQ LTSSFNPIIE VESKWYGEFY KE
Source	Yeast
Target Names	MIMI_R883
Protein Names	Recommended name: Uncharacterized protein R883
Expression Region	1-452
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.