



Recombinant Acanthamoeba polyphaga mimivirus

Uncharacterized protein R635 (MIMI_R635)

Product Code	CSB-BP733103ADAZ
Abbreviation	MIMI_R635
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UR86
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MELENICSGR VLCGNSACEK CFHKSFASHP YAKYWSENNT ESPEFVNLKS ASYFEFKCGE CHHLFSARPC YVQNKFCPYC LDVKKICVDK NCTECYHKSF ASDPLSIFWS DLNERRPRHI PKESHIRGMF NCHECGHNFE TSIAVFIKYH KIRQIICPYC SEKMLCRDKN CIKCLTASFA SRPEIVFWSS DNKMSPRHIF KNSKVNYIFS CTCGFKFTQS PLCIISSTKW WDRHDCKKYE YLKIECEKPQ YFEDFEHYLS LVSAQKSINR IQSIPYIDVF DDIQFDQQSD IQSEPSPIFT TTKSVESADD IGTSGFTESS ESTPLSDQPG CEYFEQIESV LGPDIDVDHF FSWYNETVNY NRPLKRQRQM TASEINELND PNSVYNSPEF DHQGDQKKLT EENGCVVQ
Source	Baculovirus
Target Names	MIMI_R635
Protein Names	Recommended name: Uncharacterized protein R635
Expression Region	1-408
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.