



# Recombinant *Oryza sativa* subsp. japonica Zinc finger CCCH domain-containing protein 1 (Os01g0174600, LOC\_Os01g07930)

<b>Product Code</b>	CSB-BP733149OFG
<b>Abbreviation</b>	Os01g0174600, LOC_Os01g07930
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5VR07
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAGRGGMVEW EVGRRRDESED VIVLSPGPPA RRRPPPVKAV EPESGGFAYE PPEKLFYKTR VCETFVTSGR CMFEDGCTFA HGDEELRPSL TACAGGWRKP SPSLSAAAPP VAVAPTPPPA QVVHELLARG SGSGGGGHRA ITKVCFEFRD KGICYFGETC AFPHVSAAEI RQGSRLSSMS SSSWEMPARR SVAVTVPRTF VSVPPVAPPP PPHYRVNNS SSYNAASMAA AAPAASDANL VAQQPPPEQG GRKMTRLEML SLKKMTGIYG DWLEGYEHF
<b>Source</b>	Baculovirus
<b>Target Names</b>	Os01g0174600
<b>Protein Names</b>	Recommended name: Zinc finger CCCH domain-containing protein 1 Short name= OsC3H1
<b>Expression Region</b>	1-279
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.