



# Recombinant Human Pre-mRNA-splicing factor 38B (PRPF38B)

<b>Product Code</b>	CSB-MP733154HU
<b>Abbreviation</b>	PRPF38B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5VTL8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	<p>           ANNSPALTG NSQPQHAAA AAAQQQQQCG GGGATKPAVS            GKQGNVLPW GNEKTMNLNP MILTNILSSP YFKVQLYELK TYHEVVDEIY            FKVTHVEPWE KGSRKTAGQT GMCGGVRGVG TGGIVSTAFCLLYKLFRTLK            TRKQVMGLIT HTDSPYIRAL GFMYIRYTQP PDLWDWFES FLDDEEDLDV            KAGGGCVMTI GEMLRSLTK LEWFSTLFPV IPVPVQKNID QKIKTRPRKI            KKDGGKEGAE IDRHVERRRS RSPRRSLSPR RSPRRRSRS HHREGHGSSS            FDRELEREKE RQLEREAKE REKERRRSRS IDRGLERRRS RSRRHRRSRS            RSRDRKGDRR DRDREREKEN ERGRRDRDY DKERGNEREK            ERERSRERSK EQRSRGEVEE KKHKEDKDDR RHRDDKRDSK            KEKKHSRSRS RERKHSRSR SRNAGKRSRS RSKEKSSKH NESKEKSNKR            SRSGSQGRD SVEKSKKREH SPSKEKSRKR SRSKERSHKR            DHSDSKDQSD KHDRRRSQSI EQESQEKQHK NKDETV         </p>
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PRPF38B
<b>Protein Names</b>	Recommended name: Pre-mRNA-splicing factor 38B Alternative name(s): Sarcoma antigen NY-SAR-27
<b>Expression Region</b>	2-546
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.