



Recombinant Human Pre-mRNA-splicing factor 38B (PRPF38B)

Product Code	CSB-YP733154HU
Abbreviation	PRPF38B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5VTL8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ANNSPALTG NSQPQHAAAA AAAQQQQQCG GGGATKPAVS GKQGNVLPLW GNEKTMNLNP MILTNILSSP YFKVQLYELK TYHEVVDEIY FKVTHVEPWE KGSRKTAGQT GMCGGVRGVG TGGIVSTAFCLLYKLFTLKL TRKQVMGLIT HTDSPYIRAL GFMYIRYTOP PTDLWDWFES FLDDEEDLDV KAGGGCVMTI GEMLRSLTK LEWFSTLFPR IPVVPQKNID QQIKTRPRKI KKDGGKEGAE IDRHVERRRS RSPRRSLSPR RSPRRSRRS HHREGHGSSS FDRELEREKE RQRLEREAKE REKERRRSR IDRGLERRRS RSRRHRRSR RSRDRKGDRR DRDREREKEN ERGRRDRDY DKERGNEREK ERERSRERSK EQRSRGEVEE KKHKEDKDDR RHRDDKRDSK KEKKHSRSR RERKHSRSR SRNAGKRSR RSKEKSSKH NESKEKSNKR SRSGSQGRD SVEKSKKREH SPSKEKSRKR SRSKERSHKK DHSDSKDQSD KHDRRRSQSI EQESQEKQHK NKDETV
Source	Yeast
Target Names	PRPF38B
Protein Names	Recommended name: Pre-mRNA-splicing factor 38B Alternative name(s): Sarcoma antigen NY-SAR-27
Expression Region	2-546
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.