



Recombinant Human OTU domain-containing protein 1 (OTUD1)

Product Code	CSB-EP733158HU-B
Abbreviation	OTUD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5VV17
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MQLYSSVCTH YPAGAPGPTA AAPAPPAAT PFKVSLQPPG AAGAAPEPET GECQPAAAAE HREAAVPAA KMPAFSSCFE VVSGAAAPAS AAAGPPGASC KPPLPPHYTS TAQITVRALG ADRLLLHGPD PVPGAAGSAA APRGRCLLLA PAPAAPVPPR RGSSAWLLEE LLRPDCPEPA GLDATREGPD RNFRLSEHRQ ALAAAKHRGP AATPGSPDPG PGPWGEEHLA ERGPRGWERG GDRCDAPGGD AARRPDPEAE APPAGSIEAA PSSAAEPVIV SRSDPRDEKL ALYLAEVEKQ DKYLRQRNKY RFHIIPDGNC LYRAVSKTVY GDQSLHREL EQTVHYIADH LDHFSPLIEG DVGEFIIAAA QDGAWAGYPE LLAMGQMLNV NIHLTTGGRL ESPTVSTMIH YLGPEDSLRP SIWLSWLSNG HYDAVFDHSY PNPEYDNWCK QTQVQRKRDE ELAKSMAISL SKMYIEQNAC S
Source	E.coli
Target Names	OTUD1
Protein Names	Recommended name: OTU domain-containing protein 1 Alternative name(s): DUBA-7
Expression Region	1-481
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.