



# Recombinant Human Transcription elongation regulator 1-like protein (TCERG1L)

<b>Product Code</b>	CSB-MP733161HU
<b>Abbreviation</b>	TCERG1L
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5VWI1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MQAGARFQRR RRQLQQQQPR RRQPLLWPMD AEPPIPPPWV WMVPGSAGLL RLSAGVVVPP VLLASAPPPA APLLPLPGW PAPSEPVLPL LPLPSAPDSA AAAAAHPFPA LHGQWLFGGH SPSLGLPPSS TVELVPVFPH LCPSALATPI GKSWIDKRIP NCKIFFNNSF ALDSTWIHPE ESRFFHGHEK PRLLANQVAV SLSRPAPASR PLPTVVLAPQ PIPGGCHNSL KVTSSPAIAI ATAAAAAMVS VDPENLRGPS PSSVQPRHFL TLAPIKIPLR TSPVSDTRTE RGRVARPPAL MLRAQKSRDG DKEDKEPPPM LGGGEDSTAR GNRPVASTPV PGSPWCVVWT GDDRVEFFNP TMHLSVWEKP MDLKDRGDLN RIIEDPPHKR KLEAPATDNS DGSSSEDNRE DQDVKTKRNR TEGCGSPKPE EAKREDKGTR TPPPQILLPL EERVTHFRDM LLERGVSAFS TWEKELHKIV FDPYLLLS EERKQIFEQF VKTRIKEEYK EKSKLLLLAK EEFKLLEES KVSPRTTFKE FAEKYGRDQR FRLVQKRKQD EHFFNQFILI LKKRDKENRL RLRKMR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	TCERG1L
<b>Protein Names</b>	Recommended name: Transcription elongation regulator 1-like protein
<b>Expression Region</b>	1-586
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.