



# Recombinant Buchnera aphidicola subsp. Cinara cedri 3-isopropylmalate dehydrogenase (leuB)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-BP733237BNA                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Abbreviation</b>      | leuB                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                         |
| <b>Uniprot No.</b>       | Q5WPZ9                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Immunogen Species</b> | Buchnera aphidicola subsp. Cinara cedri (strain Cc)                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Sequence</b>          | MKKKKEYRIA VLPGDGIGPE VMQEAYKILN VINKKFDFNI YTNEYDIGGI<br>AIDKYGAPLP KKTNLGCIES HAILFGSIGG KEWTHLPAD E QPERGGLLPL<br>RKFFNLFVNL RPIKLHFSLK KLSPLKKEII QDGF DILCVR ELIGGIYFGL<br>PKGKNTFDKK NISAFDTEIY STKEIEKIAD IAFQLSLKRK KKVFSIDKAN<br>VLESSMLWRE VVIQVSKRYP NVELIHLIYD NAAMQLIKNP NQFDVILCSN<br>LFGDILSDEC AALTGSLGML PSASLNESFF GLFEPAGGSA PDLKGKNIAN<br>PIAQILSLAL LFKYSLQLEK ISHLIELSVL KTLEKGYRTL DISNGKKYIT<br>TDQMGNYIAK TLSELIN |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Target Names</b>      | leuB                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Protein Names</b>     | Recommended name: 3-isopropylmalate dehydrogenase EC= 1.1.1.85<br>Alternative name(s): 3-IPM-DH Beta-IPM dehydrogenase Short name= IMDH                                                                                                                                                                                                                                                                                                    |
| <b>Expression Region</b> | 1-367                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                        |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                         |