



Recombinant Imidazolonepropionase (hutI)

Product Code	CSB-BP733345SMZ
Abbreviation	hutI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5X9L1
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MVADVLLTHF NQLFCLNDPG HPLTGQEMKK ATIVEDGYIA IKDGLIVALG SGEPDAELVG PQTIMRSYKG KIATPGI IDC HTHLVYGGSR EHEFAKKLAG VPYLDILAQG GGILSTVRAT RSASFDONLYQ KSKRLLDYML LHGVTTVEAK SGYGLDWETE KRQLDVVAAL EKDHPIDLVS TFMAAHAIPE EYKGNPKAYL DVIKDM LPV VKEKNLAEFC DIFCEKNVFT ADESRYL LSK AKEMGFKLRI HADEIASIGG VDVA AE LSAV SAEHLMMITN DGI AKLIGAG VIGNLLPATT FSLMEDTYAP ARKMIDAGMA ITLSTDSNPG SCPTANMQFV MQLGCFMLRL TPIEVLNAV T INAAYS VNRQ ERVGS LTVGK EADIAIFDAP NIDYLFYFFA TNLIHQVYKK GQLTVDRGRI L
Source	Baculovirus
Target Names	hutI
Protein Names	Recommended name: Imidazolonepropionase EC= 3.5.2.7 Alternative name(s): Imidazolone-5-propionate hydrolase
Expression Region	1-421
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.