



Recombinant Probable tRNA sulfurtransferase (thil)

Product Code	CSB-MP733370SMZ
Abbreviation	thil
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XCT0
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MDYSEIMVRH GELSTKGKNR MRFINKLKNN IQDVLAPFPA ITVRSRDRRT HVYLNQTDYQ PVVEALKLVF GVQALSPVYK LEKSVPLLVT AVQDIMTSLY RDGLTFKIAT KRSDHAFELD SRELNSLLGG AVFEVLPNIQ AQMKHPDVTL KVEIRDEAAY ISYEEIKGAG GLPVGTSKGK MLMLSGGIDS PVAGYLALKR GLDIEVVHFA SPPYTSPGAL AKAQDLTRRL TRFGGNIQFI EVPFTEIQEE IKNKAPEAYL MTLTRRFMMR ITDAIREQRK GLVIVNGESL GQVASQTLES MQAINAVTST PIIRPVVTMD KLEIHEMAQA IDTFDISIQP FEDCCTIFAP DRPKTNPGLG NAEKYEERFD IDGLVQRAVS GIVVTEITPE IVNDEVENLI DALL
Source	Mammalian cell
Target Names	thil
Protein Names	Recommended name: Probable tRNA sulfurtransferase EC= 2.8.1.4 Alternative name(s): Sulfur carrier protein ThiS sulfurtransferase Thiamine biosynthesis protein Thil tRNA 4-thiouridine synthase
Expression Region	1-404
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.