



Recombinant Ribosomal RNA small subunit methyltransferase I (rsml)

Product Code	CSB-MP733378SMZ
Abbreviation	rsml
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XDL7
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MQVQKSFKDK KTSGLTYLVP TPIGNLQDMT FRAVATLKEV DFICAEDTRN TGLLLKHFDI ATKQISFHEH NSYEKIPDLI DLLISGRSLA QVSDAGMPSI SDPGHDLVKA AIDSDIAVVA LPGASAGITA LIASGLAPQP HVFYGFLLPRK AGQQKAFFED KHHYTETQMF CESPYRIKDT LTNMLACYGD RQVVLVRELT KLFEEYQRGS ISEILSYLEE TPLKGECLLI VAGAQVDSEI ELTADVLDLV LVQKEIQAGA KPNQAIKTIA KAYQVNRQEL YQQFHDL
Source	Mammalian cell
Target Names	rsml
Protein Names	Recommended name: Ribosomal RNA small subunit methyltransferase I EC=2.1.1.198 Alternative name(s): 16S rRNA 2'-O-ribose C1402 methyltransferase rRNA (cytidine-2'-O-)-methyltransferase RsmI
Expression Region	1-287
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.